



Adjustment Guide

Use this guide for troubleshooting any concerns or difficulties you may encounter with your breast pump.



How To Readjust Pump


1. Pause the pump.
2. Separate pump from the collection kit.
3. Be sure the bottom of the diaphragm is aligned correctly - misaligning can impact suction.
4. Ensure diaphragm is back to starting position (not compressed).
5. Re-assemble pump to shield.
6. Make sure your nipple is properly aligned with the insert and shield, BEFORE pump is started. Always turn the pump off before adjusting.
7. Cover with bra, and unpause / start pump.

Loss of Suction?

It is important to properly align the nipples into the center of the shield and have the pump firmly pressed against the breast to make a suction seal. If the pump is loosely fit, the suction generated from the pump may leak and cause a low to no suction.

Leaking from Shield?

The pump needs to be at an upright angle to avoid breast milk from back flowing to the incorrect side of the shield.



Shield Sizing

Make sure that the proper shield size is being used. The nipple should not be too tight or loose inside the shield tunnel and the nipple should be able to pull back and forth freely when pumping.

It's essential that a good seals created between the breast and silicone inserts and that the correct size insert is used.