

Comparison Chart



| | | | |
|-------------------------|-----------------------------|---------------------------|----------------------------|
| BRAND | Motif | Willow | Elvie |
| MODEL | Motif Aura | Generation 3 | Pump |
| WEIGHT | .48 lbs | .76 lbs | .92 lbs |
| SUCTION STRENGTH | 280 mmHg | 245 mmHg | 220 mmHg |
| WARRANTY | 2 yr. Motor 90 day parts | 1 yr. limited warranty | 2 yr. Hubs 90 day parts |
| NOISE LEVEL | <41 dB | N/A | 31.9 dB |
| SIZING OPTIONS INCLUDED | 4 | 1 | 2 |
| USB CHARGERS INCLUDED | 2 | 1 | 1 |
| BATTERY POWER | ✓ | ✓ | ✓ |
| CLOSED SYSTEM | ✓ | ✓ | ✓ |
| AUTO-SHUTOFF | ✓ | ✓ | ✓ |
| DIGITAL DISPLAY | ✓ | ✗ | ✗ |
| BPA FREE | ✓ | ✓ | ✓ |
| DIRECT BOTTLE FEEDING | ✓ | ✗ | ✗ |
| LAST SETTING MEMORY | ✓ | ✓ | ✓ |

motif
MEDICAL